

	Frequency	Mode	Duty Cycle	Pulse Frequency	Intensity	Treatment time
Femoral Biceps Injury Acute Phase	1Mhz	Pulsed	20%	100 Hz	1.2 W/cm² (5 cm)*	1-30 min*
Rectus Femoris Injury Acute Phase	1Mhz	Pulsed	20%	100 Hz	1.0 W/cm²	1-30 min*
Anterior Tibial Injury- Acute Phase	1Mhz	Pulsed	20%	100 Hz	0.8 W/cm²	1-30 min*
Deltoid Muscle Injury- Acute Phase	1Mhz	Pulsed	20%	100 Hz	0.8 W/cm²	1-30 min*
Rhomboid Injury Acute Phase	1Mhz	Pulsed	20%	100 Hz	1.5 W/cm²	1-30 min*
Femoral Biceps Injury Subacute Phase	1Mhz	Pulsed	50%	100 Hz	1.2 W/cm²	1-30 min*
Rectus Femoris Injury Subacute Phase	1Mhz	Pulsed	50%	100 Hz	1.0 W/cm²	1-30 min*
Anterior Tibial Injury- Subacute Phase	1Mhz	Pulsed	50%	100 Hz	0.8 W/cm²	1-30 min*
Deltoid Muscle Injury- Subacute Phase	1Mhz	Pulsed	50%	100 Hz	0.8 W/cm²	1-30 min*
Rhomboid Injury Subacute Phase	1Mhz	Pulsed	50%	100 Hz	1.5 W/cm²	1-30 min*
Reduction of Muscle Spasm - Deltoid	1Mhz	Continuous	-----	-----	0.8 W/cm²	1-30 min*
Reduction of Mus- cle Spasm - Cervical Paraspinal	1Mhz	Continuous	-----	-----	0.8 W/cm²	1-30 min*
Reduction of Mus- cle Spasm - Thoracic Paraspinal	1Mhz	Continuous	-----	-----	1.5 W/cm²	1-30 min*
Reduction of Mus- cle Spasm - Lumbar Paraspinal	1Mhz	Continuous	-----	-----	1.5 W/cm²	1-30 min*
Reduction of the Muscle Spasm - Rhomboid	1Mhz	Continuous	-----	-----	1.0 W/cm²	1-30 min*
Reduction of the Muscle Spasm - Biceps or Triceps Brachii	1Mhz	Continuous	-----	-----	0.8 W/cm²	1-30 min*
Reduction of the Muscle Spasm - Gastrocnemius	1Mhz	Continuous	-----	-----	1.5 W/cm²	1-30 min*

Reduction of the Muscle Spasm - Biceps Femoris	1Mhz	Continuous	-----	-----	1.5 W/cm2	1-30 min*
Reduction of the Muscle Spasm - Rectus Femoris	1Mhz	Continuous	-----	-----	1.2 W/cm2	1-30 min*
Reduction of the Muscle Spasm - Rectus Femoris	1Mhz	Pulsed	20%	100 Hz	0.8 W/cm2	1-30 min*
Trochanteric Bursitis Chronic Phase	1Mhz	Continuous	-----	-----	0.8 W/cm2	1-30 min*
Gluteus Medius Tendinopathy Acute Phase	1Mhz	Pulsed	20%	100 Hz	1.0 W/cm2	1-30 min*
Supraspinal Tendinopathy Acute Phase	1Mhz	Pulsed	20%	100 Hz	1.0 W/cm2	1-30 min*
Gluteus Medius Tendinopathy Chronic Phase	1Mhz	Continuous	-----	-----	1.0 W/cm2	1-30 min*
Supraspinal Tendinopathy - Chronic Phase	1Mhz	Continuous	-----	-----	1.0 W/cm2	1-30 min*
Sonophoresis (Rehabilitation)	1Mhz	Continuous	-----	-----	1.0 W/cm2	1-30 min*



* Intensity calculated according to the depth of the target tissue.

** The treatment time depends of the size of the treatment area and the **ERA (Effective Radiating Area)** of the transducer. The operator must calculate the treatment time according to the size value of the treatment area divided by the value of the ERA of the transducer.