

Leading the Industry for Over 20 Years

stryker[®]
INSTRUMENTS



Power of Choice
dual trigger, dual speed
for greater versatility

System 5 Rotary Handpieces

Stryker Instruments, the leader in cordless power, brings its 5th generation of battery powered instrumentation...

SYSTEM 5

...Advancing Orthopaedics with power, speed, versatility, reliability and performance.

Collet Release

Modes of Operation:

Safe

Reverse Lockout

Forward, Reverse & Oscillate

Mode Selector Lever

4205 Dual Trigger Handpiece

- Now providing more options in rotary power handpieces.
- Offering two trigger designs for surgeon preference.



Reverse Trigger

Forward Trigger

4203 Rotary Handpiece

- Rotary handpieces allow for both high speed drilling and high torque reaming.
- 180° trigger rotation allows for quick and easy selection of direction/safety options.



Collets

Drilling & Reaming Attachments

Colored coded rings on drilling and reaming attachments allow for easy recognition and selection.


Stryker Instruments offers a wide assortment of quick connect attachments to be used with a variety of reamers, drill bits, and cutting accessories.

Drilling

Pin Collet simplifies the placement of guide pins.



4103-126 Pin Collet



4103-131 1:1 1/4" Jacob's® Drill



4103-135 Hudson® Modified Trinkle Drill



4103-160 1:1 Standard Trinkle Drill



4103-180 2:1 High Torque 1/4" Jacob's® Drill



4103-133 Keyless Drill




4103-82-131 Keyless Chuck Adapter




4103-110 1:1 Synthes® Drill


Reaming




4103-231 3.25:1 1/4" Jacob's® Reamer




4103-210 3.25:1 Synthes® Reamer




4103-213 3.25:1 Hudson® Reamer




4103-235 3.25:1 Hudson® Modified Trinkle Reamer



4103-413 4:1 High Torque Hudson® Reamer



4103-435 4:1 High Torque Hudson® /Modified Trinkle Reamer



4103-260 3.25:1 Standard Trinkle Reamer

System 5 Heavy Duty Saws


SYSTEM **5**

Advanced Motor Technology provides increased torque for a faster and more efficient cut.



4208 Sagittal Saw

- The only dual speed saw on the market.
- Two speed control allows for a wider range of cutting options.



- See bur & blade guides for a complete list of available cutting accessories.

4206 Reciprocating Saw

- New and improved motor allows for faster cutting action.



4207 Reciprocating/Sternum Saw


- Extended run battery available for higher speed and longer run time.

Accessories

4102-450 Sterilization Case

- Holds 3 handpieces and all accessories for complete sterilization.

- Convenient inner tray presents instruments on a sterile field.




4110-120 Battery Charger

- Quickly charges System 5 and Cordless Driver 2 batteries.
- Non-sterile battery (red) available.
- Conditioning cycle completely discharges non-sterile battery before charging.
- Modular design allows for 4 different batteries to be charged at the same time.
- Indicator lights provide battery status feedback with "Charge", "Ready", and "Replace" indicators.



- Battery charging station diagnoses battery levels and usage.



Complete System 5 Ordering Information

System 5 Handpieces

PART #	DESCRIPTION
4203	Rotary Handpiece, Single Trigger
4205	Rotary Handpiece, Dual Trigger
4206	Reciprocating Saw
4208	Sagittal Saw

Drilling Attachments

PART #	DESCRIPTION
4103-82-131	Keyless Jacobs® Chuck Adapter
4103-110	Synthes® Chuck
4103-126	Pin Collet
4103-131	1/4" Drill with Jacobs® Chuck (1:1)
4103-133	Keyless Jacobs® Drill Attachment (1:1)
4103-135	Hudson®/Modified Trinkle Drill
4103-160	Trinkle Drill Attachment (1:1)
4103-180	High Torque 1/4" Drill with Jacobs® Chuck (2:1)

Reaming Attachments


PART #	DESCRIPTION
4103-210	Synthes® Reamer (3.25:1)
4103-213	Hudson® Reamer (3.25:1)
4103-231	1/4" Reamer with Jacobs® Chuck (3.25:1)
4103-235	Hudson®/Modified Trinkle Reamer (3.25:1)
4103-260	Standard Trinkle Reamer (3.25:1)
4103-413	High Torque Hudson® Adapter
4103-435	High Torque Hudson®/Modified Trinkle Reamer (4:1)

System 5 Accessories

PART #	DESCRIPTION
4115	System 5 Battery
4126	System 5 Aseptic Battery Kit (includes battery pack, housing, & shield)
4126-110	Non-Sterile Large Battery
4126-120	System 5 Aseptic Battery Housing
4126-130	Large Transfer Shield 4110-115 Battery Adapter (for use with 4115 batteries in System 2000 charger)
4110-120	Battery Charger, 120 Volt
4110-215	Charger Module (2115 battery in System 5 Charger)
4110-415	Charger Module (4115 battery)
4110-426	Non-Sterile Battery Charger Module (4126-110 battery)
5100-901-906	Electrical Connection Brush
2106-151	Sterilization Case (holds 4 handpieces)
2106-152	Sterilization Basket (holds 2 handpieces)
4102-450	Plastic Sterilization Case (holds 4 handpieces)
4102-451	Sterilization Case (holds 3 handpieces)
4102-452	Mini Sterilization Case (holds 2 handpieces)

Sternum Applications

PART #	DESCRIPTION
4207	Reciprocating/Sternum Saw
4112	Cordless Driver Battery (for use with 4207 Sternum Saw)
4212	Small Smart Cordless Driver Battery (for use with 4207)
4215	Large Smart Cordless Driver Battery (for use with 4207)
4110-412	Charger Module (4112 battery)
4110-422	Cordless Driver Smart Charger Module (4212 & 4215 batteries)
4107-008	Sternum Blade Guard
298-97-100	Sternum Blade*
298-97-101	Aggressive Sternum Blade



4126 System 5 Aseptic Battery Kit
(includes battery pack, housing, & shield)

- Non-sterile battery ensures longer life of battery.
- Battery is easily charged in System 5 charger.
- Simple assembly for operating room staff.



Battery Pack



Housing



Shield

* See Bur and Blade literature for complete listing of available cutting accessories