

BONE MARROW BIOPSY NEEDLES



BONE MARROW BIOPSY NEEDLES

The "Double Diamond" Stylet Point is designed for quick and easy penetration into the marrow cavity.

- Ergonomically designed handle allows the clinician to apply sufficient pressure for penetration without excessive hand fatigue.
- Luer Lock Connector on the handle provides for secure syringe attachment.
- The Marrow Extraction Cannula included with our BMNJ-MEC-1110 provides greater patient comfort, easier marrow sampling, low risk of sample loss and better quality specimens.

BMNJ-MEC-1110	11g x 10cm with Extraction Cannula	Box of 10
JD-1106	11g x 15cm	Box of 10

BONE MARROW ASPIRATION NEEDLES

The combination of an ergonomically designed handle with an exceptionally sharp needle point provides the clinician with comfort, while easily penetrating hard bones.

- Adjustable needle stop ensures safe aspirations at the sternum by controlling the depth of penetration.
- The needle stop can be removed if aspirating at the iliac crest.

1501-I	15g x 5cm	Box of 10
1601-I	16g x 5cm	Box of 10

**Volume discounts and price matching available.
Free samples provided upon request.**

1-800-620-2101 www.mpmmedicalsupply.com



265 Willow Brook Road Unit 10, Freehold, NJ 07728 Phone (800) 620-2101 Fax: (800) 620-2127
Email: Sales@mpmmmedicalsupply.com

HEMATOLOGY ONCOLOGY



**Your Premium Source for
Quality Hematology Oncology
Products & Instruments**



**Visit Us Online!
www.mpmmedicalsupply.com**

BONE MARROW BIOPSY TRAYS

BONE MARROW BIOPSY TRAYS

BONE MARROW BIOPSY TRAYS

Bone marrow biopsy trays are designed to provide the majority of components necessary to perform bone marrow biopsy procedures.



BMTX-1023 11g x 10cm "J" needle and 15g x 5cm "I" needle Box of 10

TRAY INCLUDES:

- | | |
|-----------------------------------|--------------------------------------|
| (1) BMB "J" Needle, 11g x 10cm | (1) Lidocaine Insert |
| (1) BMB "I" Needle, 15g x 5cm | (10) Glass Slides, Frosted |
| (1) Obturator | (1) Drape, Fenestrated with Adhesive |
| (1) Cap – Luer Lock, Red | (10) Guaze, 3"x3" |
| (1) Anti Stick Probe Guide, Green | (1) Adhesive Bandage |
| (1) Filter Needle 19g x 1-1/2" | (1) Towel |
| (1) Needle, 21g x 1-1/2" | (1) Scalpel, #11 Mini |
| (1) Needle, 25g x 5/8" | (3) PVP Swabsticks |
| (1) Syringe, 5cc L/L | (1) Label, Patient |
| (1) Syringe, 20cc L/L | (1) CSR Wrap 30" |
| (2) Lidocaine 1% 5cc | |



BMTX-1024-MEC 11g x 10cm "J" needle with extraction cannula and 15g x 5cm "I" needle Box of 10

TRAY INCLUDES:

- | | |
|-----------------------------------|--------------------------------------|
| (1) BMB "J" Needle, 11g x 10cm | (2) Lidocaine 1% 5cc |
| (1) BMB "I" Needle, 15g x 5cm | (1) Lidocaine Insert |
| (1) Marrow Extraction Cannula | (10) Glass Slides, Frosted |
| (1) Obturator | (1) Drape, Fenestrated with Adhesive |
| (1) Cap – Luer Lock, Red | (10) Guaze, 3"x3" |
| (1) Anti Stick Probe Guide, Green | (1) Adhesive Bandage |
| (1) Filter Needle 19g x 1-1/2" | (1) Towel |
| (1) Needle, 21g x 1-1/2" | (1) Scalpel, #11 Mini |
| (1) Needle, 25g x 5/8" | (3) PVP Swabsticks |
| (1) Syringe, 5cc L/L | (1) Label, Patient |
| (1) Syringe, 20cc L/L | (1) CSR Wrap 30" |



BMTX-1027 11g x 10cm "J" needle Box of 10

TRAY INCLUDES:

- | | |
|-----------------------------------|--------------------------------------|
| (1) BMB "J" Needle, 11g x 10cm | (1) Lidocaine 1% 5cc |
| (1) Obturator | (1) Lidocaine Insert |
| (1) Cap – Luer Lock, Red | (1) Drape, Fenestrated with Adhesive |
| (1) Anti Stick Probe Guide, Green | (10) Guaze, 4"x4" |
| (1) Needle, 21g x 1-1/2" | (1) Scalpel, #11 Mini |
| (1) Needle, 20g x 1-1/2" | (3) PVP Swabsticks |
| (1) Needle, 25g x 5/8" | (1) Foam Ampule Holder |
| (1) Syringe, 10cc L/L | (1) CSR Wrap 24" |
| (1) Syringe, 20cc L/L | |



BMTX-1028 11g x 10cm "J" needle and 15g x 5cm "I" needle - Economy Box of 10

TRAY INCLUDES:

- | | |
|-----------------------------------|--------------------------------------|
| (1) BMB "J" Needle, 11g x 10cm | (1) Syringe, 5cc L/L |
| (1) BMB "I" Needle, 15g x 5cm | (1) Syringe, 20cc L/L |
| (1) Obturator | (1) Drape, Fenestrated with Adhesive |
| (1) Cap – Luer Lock, Red | (5) Guaze, 3"x3" |
| (1) Anti Stick Probe Guide, Green | (1) Scalpel, #11 Mini |
| (1) Needle, 21g x 1-1/2" | (3) PVP Swabsticks |
| (1) Needle, 20g x 1-1/2" | (1) Drape/Towel, 19" x 26" 3-ply |
| (1) Needle, 25g x 5/8" | (1) CSR Wrap 24" |



BMT1110J-151I-MS 11g x 10cm "J" needle w/ extraction cannula and 15g x 5cm "I" needle with Sharps Safe components Box of 10

TRAY INCLUDES:

- | | |
|-----------------------------------|--------------------------------------|
| (1) BMB "J" Needle, 11g x 10cm | (1) Lidocaine Insert |
| (1) BMB "I" Needle, 15g x 5cm | (10) Glass Slides, Frosted |
| (1) Marrow Extraction Cannula | (1) Drape, Fenestrated with Adhesive |
| (1) Obturator | (10) Guaze, 3"x3" |
| (1) Cap – Luer Lock, Red | (1) Adhesive Bandage |
| (1) Anti Stick Probe Guide, Green | (1) Towel |
| (1) Filter Needle 19g x 1-1/2" | (1) Safety Scalpel, #11 |
| (1) Safety Needle, 21g x 1-1/2" | (3) Point-Lok Safety Device |
| (1) Safety Needle, 25g x 1-1/2" | (3) PVP Swabsticks |
| (1) Syringe, 10cc L/L | (1) Label, Patient |
| (1) Syringe, 20cc L/L | (1) CSR Wrap 30" |
| (2) Lidocaine 1% 5cc | |

CUSTOMIZE YOUR TRAY

Trays can be customized to meet the specifications of your facility by varying the quantity of standard components, eliminating certain items or adding new componentry.