

8.1200°C -1700°C Dual Heating Zone Tube Furnace



Model	Max. temp	Operating Max temp.	Heating length	Chamber size (Φ x Length)	Voltage	Power	Heating element	Temp. Precision
STGS-40-12	1200°C	1150°C	600mm	Φ40×120mm	380 V	6KW	Electric Alloy Wire	±1°C
STGS-60-12			600mm	Φ60×120mm	380 V	6KW		
STGS-80-12			600mm	Φ80×120mm	380 V	8KW		
STGS-100-12			600mm	Φ100×120mm	380V	10KW		
STGS-40-14	1400°C	1350°C	600mm	Φ40×120mm	380 V	6KW	SiC Rod	±1°C
STGS-60-14			600mm	Φ60×120mm	380 V	6KW		
STGS-80-14			600mm	Φ80×120mm	380 V	8KW		
STGS-100-14			600mm	Φ100×120mm	380 V	10KW		
STGS-60-17	1700°C	1650°C	600mm	Φ60×120mm	380 V	12KW	MoSi2 Rod	±1°C