

Tissue Floating Bath

Model No.: XH-1001



Feature

- This instrument has been developed using microprocessor control and Modular . Design that has improved its reliability and precision, and it will be reliable for your lab.
- Material: steel and with white printing coating. The anodized black surface offers best contrast in identifying sections. which makes selecting and picking up sections onto a slide efficient and convenient. The transparent glass lid keeps the heat more evenly and reduces heat loss.

Technical Specification

Outside Dimensions	ø270mm, H210mm
Bath Dimensions	Ø 210 mm, Depth 65 mm
Bath temperature range	ambient - 75°C
Capacity	4L
Voltage	110V, 60Hz / 220V, 50Hz
Power	400W
Relative humidity	Less then 85%