

9. Vacuum Split Tube Furnace





Model	Max. temp	Operating Max temp.	Heating Length	Tube size (Φ x Length)	Voltage	Power	Tube material	Heating Element
STGK-40-12	1200°C	1150°C	300 mm	Φ40×100mm	220V	2KW	Quartz Tube	±1°C
STGK-60-12			300 mm	Φ60×100mm	220V	3KW		
STGK-80-12			300 mm	Φ80×100mm	220V	4KW		
STGK-100-12			300 mm	Φ100×100mm	220V	5KW		