


Transparent Water Bath

◎ Used for precision constant temperature and auxiliary.



Characteristics:

- Using innovative design structure, effectively prevent steam leakage.
- Intelligent temperature control instrument with high brightness digital display, light-touch type key has the advantages of simple operation and high precision temperature control.
- Shell, interior and cover of the product are all made of transparent engineering plastics, can clearly see the object to be heated.

Main Technical Parameters:

Model	DK-98-IV
Voltage (V)	200-240V/100-120V
Power (kW)	0.8
Temp. Motion (°C)	± 0.5
Temp. Range (°C)	Rt+5~100
Temp. Sensitivity (°C)	≤ ± 1
Shelves Height (mm)	Distance from bottom 100mm
Shelves Holes Diameter	Ø16-48mm
Inner Chamber Size (L*W*H)	280 x 220 x150mm
Exterior Size (L*W*H)	400 x260x 190mm
Packing Size (L*W*H)	450 x350 x270mm
N.W. (KG)	3.5
G.W. (KG)	5
Optional accessories	Different size chamber