

Vacuum Drying Oven



Characteristics:

- 4-side heating method, temperature conduction by shelves, more efficiency.
- With inflatable port, once-forming thermostability silica gel sealant strip, with good tightness
- Double observation glass window, inner is explosion-proof glass, outside is resin protection plate, more safely.
- Fastening handle, convenient for opening and closing.

Main technical parameters:

Model		DZ-1A II	DZ-2A II	DZ-3A II	
		DZ-1BC II	DZ-2BC II	DZ-3BC II	
Heating mode		Four walls heat conduction and decompression			
Function	Temp.Range	RT+10-250°C			
	Useful vacuum degree Range	<133Pa			
	Temp.Resolution Ratio	0.1°C			
	Temp.motion	±1°C			
	Heating-up time	80min	100min	120min	
Structure	Inner Chamber	BCmodel:High quality hairline stainless steel sheet;A model:Anticorrosion treatment of high strength galvanized sheet			
	Outer Shell	Cold rolling steel electrostatic spraying exterior			
	Insulation layer	Aluminum silicate fiber			
	Heater	Stainless steel electric heating tube			
	Observation window	Bullet proof glass,acrylic exterior protection			
	Vacuum meter	Pointer type,accuracy calss 2.5			
	Nozzle diameter	10mm			
	Power rating	0.8kW	1.4kW	2.0kW	
	Controller	Temp.control mode	Two temoerature section PID intelligent		
		Temp.setting mode	Touch button setting		
Temp.display mode		Measuring temperature:display on the four digital tubes upper row;setting temperature:display on the lower row			
Timer		0-9999min(with timing wait function)			
Operation function		Fixed temperature operation,timing function,auto stop			
Additional function		Sensor deviation correction,temperature overshoot self-tuning,internal parameter locking,power-off parameter memory			
Sensor		pt100			
Safety device		Over temoerature sound-light alarm			
Specification	Inner Chamber size(W*L*H) (mm)	300*300*270	415*370*340	450*450*450	
	Exterior size(W*L*H)(mm)	480*480*606	560*540*680	640*590*780	
	Packing size(W*L*H)(mm)	590*550*750	704*620*814	739*700*889	
	Volume	24L	52L	91L	
	Number of shelf	2			
	Load per rack	15kg			
	Shelf space	100mm	140mm	185mm	
	(50/60HZ)Current rating	AC220V/3.6A	AC220V/6.3A	AC220V/9.1A	
	NW/GW (kg)	42/52	67/92	80/105	
Accessory	Shelf	2			
	Shelf frame	4			
Optional accessories		Silica gel connecting pipe,vacuum pump			