


<b>ITEM:</b>	<b>DC Series Water Bath</b>	
<b>PART NUMBER:</b>	DC-4006	
<b>BASIC FUNCTION:</b>		
<p>The series of Refrigerated thermostatic bath has adopted the latest international advanced CFC-free refrigeration system. The main components used are imported and the performance is stable and reliable. It is widely used in research institutes, universities, corporate quality control departments and industrial sectors relating to petroleum, chemical engineering, electronic instruments, physics, chemistry, bio-engineering, medicine and health, life science, light industry and food, materials testing and chemical analysis etc. It is used to put trial samples or newly-made products to a constant-temperature test. It can also be used as a heat source or cold source for a direct heating or cooling and auxiliary heating or cooling.</p>		
<b>SPECIFICATIONS: (As per one set)</b>		
Temperature Range (°C):	-40~100	
Temperature Stability (°C) :	±0.1	
Display Resolution (mm3):	0.1/0.01	
Inner Tank size (mm3)	260×200×140	
Depth (mm)	140	
Pump Flow (L/min)	6	
Top Opening Dimen (mm)	180×140	
Drain	√	
N.W (kg)	36	
Packing Size (mm)	510*480*1015	