



## HYDROMASSAGE TABLE

### GOOD FOR YOUR PATIENTS

Enhance your patient's physical and mental well-being naturally, using hydromassage. The pressure of the heated water impacting the body releases tension from the muscles. As the body relaxes, stress is relieved, and the mind relaxes.

### GREAT FOR YOUR PRACTICE

Hydromassage improves patient satisfaction. Once they try it, they will ask for it every time. Happy, satisfied patients keep their appointments, reschedule more often, and refer their family and friends.


Manufactured By: Sidmar, PO Box 249, 31530 125th Street NW, Princeton, MN 55371 USA

# What is Hydromassage Therapy?

Imagine how you felt the last time you were in a hot tub. The warmth of the water relaxed you and the pressure of the water jets against your body soothed your aches and pains. Now imagine feeling this way without undressing or getting wet. That's what our hydromassage tables feel like. Just lie back, set the timer and within fifteen minutes you will achieve a state of total relaxation.

## How it Helps Your Body

Warm water massage increases blood flow, which speeds the delivery of oxygen and vital nutrients, aiding and speeding recovery. As your muscles warm, they relax. Muscle pain is reduced as increased blood flow carries away toxins in tired, aching muscle tissue. The temperature increase stimulates the release of endorphins, which work to naturally relieve pain, and improve your sense of well-being.



## Benefits of Hydromassage Therapy for Your Patients

- Diminishes aches and pains
- Relaxed muscles move easier with increased range of motion
- Increased blood circulation improves body function
- Overall relaxation and calmness, reduces stress and clears the mind
- Relaxed patients respond better to your care

## Specifications

<b>Model:</b>	Pro S10
<b>Upholstery:</b>	Naugahyde Spirit Millennium Vinyl - Attractive, strong, and durable.
<b>Color Choices:</b>	black, burgundy, brown, navy, tan, teal, green, gray, blueridge
<b>Controls:</b>	Session timer (3 to 30 minutes), temperature (70 to 100 F), and massage intensity (12 psi to 17 psi)
<b>Dimensions:</b>	73" L x 23" W x 23" H
<b>Water Volume:</b>	21 gallons
<b>Dry Weight:</b>	132 lbs.
<b>Filled Weight:</b>	315 lbs.
<b>Electrical:</b>	120 volts, 14 amps